

Defiant Technologies' Bubble Flow Meter

Know Your Flow



**Measures flow rates from
pressure and vacuum**


Features

- ✓ Scale markings at 0.5 mL, 1 mL, 5 mL, 10 mL, 25 mL, and 50 mL
- ✓ Convenient clamp stand for mounting on bench, desk, or table
- ✓ Adapter for measuring flow into vacuum

Accurate Flow Measurement Made Easy

Accurate flow measurements are essential for reliable chromatography.

- ✓ During sample collection, flow rates determine how much of analyte is collected for analysis
- ✓ During analysis, flow rates govern the time analytes emerge from the column
- ✓ Measure the flow coming and going with the three stage bubble flow meter from Defiant Technologies



Includes:

- Three-stage glass tube for wide measurement range
- Hardened case for hauling into the field
- Clamp stand for quick set-up anywhere
- Top adapter for measuring suction flows
- Stop watch not included (use the one on your phone, it's OK)

Defiant Technologies, Inc.

6814A Academy Parkway W NE

Albuquerque, NM 87109

www.defiant-tech.com

(505) 999-5880

sales@defiant-tech.com

Publication 20180814